

Product Datasheet

Canine Growth Hormone R (GHR) Protein, His Tag (orb1818956)

Description

Canine Growth Hormone R, His Tag is expressed

Reactivity	Canine
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Growth Hormone R
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 29.2 kDa. The protein migrates as 43-53 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Canine Growth Hormone R, His Tag (orb1818956) is expressed from human 293 cells (HEK293). It contains AA Thr 27 - Gln 264 (Accession # Q9TU69-1).
MW	29.2 kDa
Uniprot ID	Q9TU69
Expiration Date	6 months from date of receipt.



Loaded
Canine
Growth
Hormone
R, His Tag
...