
Product Datasheet

**Human Erythropoietin R / EPO R Protein, Fc Tag
(MALS verified) (orb1818958)**

Description

Human Erythropoietin R Protein,

293 cells (HEK293). It contains AA Al...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Erythropoietin R

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 51.3 kDa. The protein migrates as 58-63 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Human Erythropoietin R Protein, Fc Tag (orb1818958) is expressed from human 293 cells (HEK293). It contains AA Ala 25 - Pro 250 (Accession # P19235-1).

MW

51.3 kDa

Uniprot ID
P19235
Expiration Date

6 months from date of receipt.