
Product Datasheet

Anti-XDH Antibody (orb1819408)

Description

Anti-XDH Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, FC, IHC, WB

Immunogen

E.coli-derived human XDH recombinant protein (Position: R439-R794). Human XDH shares 90.4% amino acid (aa) sequence identity with mouse XDH.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Western blot, 0.25-0.5 µg/ml, Mouse, Rat Immunohistochemistry, 2-5 µg/ml, Human Flow Cytometry (Fixed), 1-3 µg/1x10⁶ cells, Human ELISA, 0.1-0.5 µg/ml, Human. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml

Clonality

Polyclonal

MW

150 kDa

Uniprot ID
P47989
Expiration Date

12 months from date of receipt.