

Product Datasheet

Mouse IL-36 gamma Recombinant Protein (orb1819504)

Description	The Mouse IL-36 gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-36 gamma
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GRETPDFGEV FDLDQQVWIF RNQALVTVPR SHRVTGPSVT ILPCKYPESL EQDKGIAIYL GIQNPDKCLF CKEVNGHPTL LLKEEKILD YHHPEPMKPF LFYHTRTGGT STFESVAFPG HYIASSKTGN PIFLTSSKGE YYNINFNLDI K (151)
Purity	98%
Source	Yeast
MW	17.2 kDa
Entrez	215257
Expiration Date	6 months from date of receipt.