

Product Datasheet

Ferret IL-5 Recombinant Protein (orb1819509)

Description	The Ferret IL-5 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-5
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	TPVEDPMNRL VAETLTLLST HRTLLIGDGN LMIPTPENKN HQLCIEEVFQ GIDILKNRTA QAAAVDKLFQ NLSFIKAHID SQKKKCGGER WRVKKFLDYL QVFLGVINTE WTES (115)
Purity	98%
Source	Yeast
MW	13.1 kDa
Entrez	101688037
Expiration Date	6 months from date of receipt.