

Product Datasheet

Canine CCL8 Recombinant Protein (orb1819514)

Description	The Canine CCL8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QPDSVSIPIT CCFSMVKRKI PMQKLESYMR ITNSQCPQEA VIFKTKASRE ICADPKQKWV QDYMNHLDQK SQAQKP (76)
Purity	98%
Source	Yeast
MW	8.8 kDa
Entrez	448792
Expiration Date	6 months from date of receipt.