

Product Datasheet

Mouse IL-36RA Recombinant Protein (orb1819521)

Description	The Mouse IL-36 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturally...
Target	IL-36RA
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VLSGALCFRM KDSALKVLYL HNNQLLAGGL HAEKVIK GEE ISVVPNRALD ASLSPVILGV QGGSQLSCG TEKGPIKLE PVNIMELYLG AKESKSFTFY RRDMGLTSSF ESAAYPGWFL CTSPEADQPV RLTQIPEDPA WDAPITDFYF QQCD (154)
Purity	98%
Source	Yeast
MW	16.9 kDa
Entrez	54450
Expiration Date	6 months from date of receipt.