

Product Datasheet

Human IL-36RA Recombinant Protein (orb1819526)

Description	The Human IL-36 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturally...
Target	IL-36RA
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MVLSGALCFR MKDSALKVLY LHNNQLLAGG LHAGKVIKGE EISVVPNRWL DASLSPVILG VQGGSQLCSC GVGQEPTLTL EPVNIMELYL GAKESKSFTF YRRDMGLTSS FESAAYPGWF LCTVPEADQP VRLTQLPENG GWNAPITDFY FQQC (154)
Purity	98%
Source	Yeast
MW	16.8 kDa
Entrez	26525
Expiration Date	6 months from date of receipt.