

Product Datasheet

Canine IL-36RA Recombinant Protein (orb1819529)

| | |
|-------------------------|--|
| Description | The Canine IL-36 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturall... |
| Target | IL-36RA |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | VLSGALCFRM KDAALKVLYL HDNQLLAGGL HAGKVIKGE ISVVPNRSLD AKLSPVILGV QGGSQLSCG TGQEPTLKE PVNIMELYLG AKESKSFTFY RRDTGLTSSF ESAAYPGWFL CTVPEADQPL RLTQLSEDAS WDNPITDFYF QQCD (154) |
| Purity | 98% |
| Source | Yeast |
| MW | 16.9 kDa |
| Entrez | 611869 |
| Expiration Date | 6 months from date of receipt. |