

## Product Datasheet

### Feline IL-36RA Recombinant Protein (orb1819530)

<b>Description</b>	The Feline IL-36 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturall...
<b>Target</b>	IL-36RA
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	VLSGALCFRM KDAELKVLVYL HDNQLLAGRL HAGNIIQGEE ISVVPNRS LD AKLSPVILGV QGGSQCLSCG TGQEPTLKLE PVNIMELYLG AKESKSFTFY RRDTGLTSSF ESAAYPGWFL CTVPEADQPL RLTQLSEDAS WDSPITDFYF QQCD (154)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	16.9 kDa
<b>Entrez</b>	<b>101097230</b>
<b>Expiration Date</b>	6 months from date of receipt.