

Product Datasheet

Feline IL-36RA Recombinant Protein (orb1819532)

Description	The Feline IL-36 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturall...
Target	IL-36RA
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VLSGALCFRM KDAELKVLVYL HDNQLLAGRL HAGNIIQGEE ISVVPNRS LD AKLSPVILGV QGGSQCLSCG TGQEPTLKLE PVNIMELYLG AKESKSFTFY RRDTGLTSSF ESAAYPGWFL CTVPEADQPL RLTQLSE DAS WDSPITDFYF QQCD (154)
Purity	98%
Source	Yeast
MW	16.9 kDa
Entrez	101097230
Expiration Date	6 months from date of receipt.