

Product Datasheet

Equine CXCL3 Recombinant Protein (orb1819535)

Description	The Equine CXCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CXCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APVVSELRCQ CLQTVQGIHQ KNIQSVKVTP AGSHCAQTEV IATLKNGQKA CLNPAAPMVK KMIDRMLNKG STN (73)
Purity	98%
Source	Yeast
MW	7.8 kDa
Entrez	100056258
Expiration Date	6 months from date of receipt.