

Product Datasheet

Feline CCL17 Recombinant Protein (orb1819539)

Description	The Feline CCL17 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CCL17
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ARATNVGRECCLEYFKGAIPLRRLTGWYRTSGECSKDAIVFVTIHGKSICSDPKDT RVKKTVRYLQSIMNPVPE (75)
Purity	98%
Source	Yeast
MW	8.5 kDa
Entrez	493832
Expiration Date	6 months from date of receipt.