

## Product Datasheet

### Guinea Pig IL-6 Recombinant Protein (orb1819540)

<b>Description</b>	The Guinea Pig IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
<b>Target</b>	IL-6
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	AFPTAQVQHD FTADTTDEMT TAEMTTTTPN KPTTSASQVF QMFMRVYQAV KELKNEMNKH NVEKAILNNL DLPKLELDG CFFNGYNWET CQLKITPGLF KFQTYLQSMQ NKLQNESENK KAANIYAGIK SLSLFMFSKI NNTEQMEFLS PTPDATLLEK LETQSQTQML LIAEIVLQRL EEFLQDSLRA IRKADWEGRN (200)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	23.0 kDa
<b>Entrez</b>	<a href="#">106027696</a>
<b>Expiration Date</b>	6 months from date of receipt.