

Product Datasheet

Bat FLT3 Ligand Recombinant Protein (orb1819541)

Description	The Jamaican fruit-eating Bat FLT3 Ligand yeast-derived recombinant protein is not tagged and is nat...
Target	FLT3 Ligand
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	RTPDCTFPHS PISSTFADTI RKLSDYLLD YPVTVATNLH DDKLCWAFWR LVLAQRCMER LKTVAGSQMQ QLLEDVNTEI YFVTSCALQD TSEQLVALKP WITRWNFSGC LELRCQPDSS TLLPSRSPKA LEATALPAAQ AP (142)
Purity	98%
Source	Yeast
MW	16.0 kDa
Entrez	119046007
Expiration Date	6 months from date of receipt.