

Product Datasheet

Rabbit TWEAK (TNFSF12) Recombinant Protein (orb1819547)

Description	The Rabbit TWEAK yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	TWEAK
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AAHYEVHPRP GQDGAQAGVD GTVSGWEEAK INSSNPLRYD RQIGKFIVTR AGLYYLYCQV HFDEGKAVYL KLDLLVDEVL ALRCLEEFSA TAASSLGPQL RLCQVSGLLP LRPGSSLRIR TLPWAHLKAA PFLTYFGLFQ VH (142)
Purity	98%
Source	Yeast
MW	15.7 kDa
Entrez	100356651
Expiration Date	6 months from date of receipt.