

Product Datasheet

Bovine IL-36 alpha Recombinant Protein (orb1819559)

Description	The Bovine IL-36 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
Target	IL-36 alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	NIQDINHLVW VLQGQALTAV PRKDQMDPVT VTLVSCKYTE TLEKGRGNPV YLGLKEPELC LFCTKVKGQP TLQLQERNIM DLYHQTEPVK PFLFYHDQNG RASSFESVAF PGWFIGSCSC GGCPVIITQE LGKIYTTDFG FTVLQP (146)
Purity	98%
Source	Yeast
MW	16.4 kDa
Entrez	523429
Expiration Date	6 months from date of receipt.