

Product Datasheet

Hamster CCL5 Recombinant Protein (orb1819569)

Description	The Hamster CCL5 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CCL5
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPRGSDDTPC CFDYLSFTLP RVHIKEYFYT SSKCSNFAVV FVTRRRNRQVC ANPKKKWVQE YINYLEMK (68)
Purity	98%
Source	Yeast
MW	8.1 kDa
Entrez	101832468
Expiration Date	6 months from date of receipt.