

Product Datasheet

Chicken IL-3 Recombinant Protein (orb1819570)

Description	The Chicken IL-3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPLDKDCP STSSLSEVIE DVEKLQGSNE AIPPNNIILM MEQGEENAVT YLEEFIRALE SSKKQIKKDF LRNLQKIYST GKNCTGWLIH NWEPMESDDF FKNFEELLKF LNMYGVSKKT S (121)
Purity	98%
Source	Yeast
MW	13.9 kDa
Entrez	474356
Expiration Date	6 months from date of receipt.