

Product Datasheet

Rabbit IL-36 alpha Recombinant Protein (orb1819574)

Description	The Rabbit IL-36 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
Target	IL-36 alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	STPAPVLRHI QDLNHQVWVW QDQTLTAVPR KHSMPATVS LIPCQYPETL ETEKGIPMYL GLDKPDVCLF CTKTREQPML QLEDKNIMTL YHHPEPVKPF LFYHSGSRT STFESVAFPG WFIAVCSEGG CPLFLTQDLG KAYTTDFELT PVN (153)
Purity	98%
Source	Yeast
MW	18.3 kDa
Entrez	100355520
Expiration Date	6 months from date of receipt.