

Product Datasheet

Hamster CCL3 Recombinant Protein (orb1819578)

Description	The Hamster CCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APYGADTLTS CCFFYGRQIS RKFIADYFET SSLCSEPGVI FLTQRNRQFC ADPKETWVQE YVADLELNAG EASSQ (75)
Purity	98%
Source	Yeast
MW	8.5 kDa
Entrez	101832203
Expiration Date	6 months from date of receipt.