

Product Datasheet

Cynomolgus Monkey IL-16 Recombinant Protein (orb1819579)

Description	The Cynomolgus Monkey IL-16 yeast-derived recombinant protein is not tagged and is naturally endotox...
Target	IL-16
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	PDLNSTTDSA ASASAASDVS VESSAEATVY TVTLEKMSAG LGFSLEGGKG SLHGDKPLTI NRIFKGAASE QSETIQPGDE ILQLAGTAMQ GLTRFEAWNI IKALPDGPVT TVIRRKSLQP KETTAAADS (129)
Purity	98%
Source	Yeast
MW	13.4 kDa
Entrez	574100
Expiration Date	6 months from date of receipt.