

Product Datasheet

Rat IL-6 Recombinant Protein (orb1819592)

Description	The Rat IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The pu...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FPTSQVRRGD FTEDTTHNRP VYTTSQVGGL ITYVLREILE MRKELCNGNS DCMNSDDALS ENNLKLPEIQ RNDGCFQGY NQEICLLKIC SGLLEFRFYL EFVKNNLQDN KKDKARVIQS NTETLVHIFK QEIKDSYKIV LPTPTSNALL MEKLESQKEW LRTKTIQLIL KALEEFLKVT MRSTRQT (187)
Purity	98%
Source	Yeast
MW	21.7 kDa
Entrez	24498
Expiration Date	6 months from date of receipt.