

Product Datasheet

Dolphin CCL24 Recombinant Protein (orb1819595)

Description	The Dolphin CCL24 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	CCL24
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MIIPSSCCLA FISKRIPERR VINYQLSNRS VCPKAGVIFT TRKGQRFCGN PKLLWVQRYM KNLDAKQKKA SSRTRVMGAK VAVQRHPANS TFI (93)
Purity	98%
Source	Yeast
MW	10.6 kDa
Entrez	101331501
Expiration Date	6 months from date of receipt.