

Product Datasheet

Dolphin CCL17 Recombinant Protein (orb1819605)

Description	The Dolphin CCL17 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	CCL17
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ARATNVGWEC CWEYFKGAIP FRKLVRWYRT PEDCPRDAIV LVT LHGRSIC SDPKDRWVKK AVKYLQNTTK SCDLAAQPA (79)
Purity	98%
Source	Yeast
MW	9.1 kDa
Entrez	101325821
Expiration Date	6 months from date of receipt.