

Product Datasheet

Rabbit CTLA-4 Recombinant Protein (orb1819615)

Description	The Rabbit sCTLA-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	CTLA-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LHVSQPAVVLASSRGVASFVCEYASSHKATEVRVTVLRQANSQMTEVCAMTYTV ENELTFIDDSTCTGISHGNKVNLTIQGLSAMDTGLYICKVELMYPPPPYYVGMGNG TQIYVIEPEPCPDSD (124)
Purity	98%
Source	Yeast
MW	13.5 kDa
Entrez	100009412
Expiration Date	6 months from date of receipt.