

Product Datasheet

Ferret IL-22 Recombinant Protein (orb1819624)

Description	The Ferret IL-22 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-22
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	EAVPIGSHCR LANSDFEQPY ITNRTFMLAG EASLADNNTD VRLIGDTLFH GVNKEDHCYL MKEVLNFTLE EVLLPQSHRF QPVMQEVESF LTKLSNRLSQ CHIDSDDQHI QRNVQNLKDT VKKLGEDGEI KAIGELNLLF IALRDACI (148)
Purity	98%
Source	Yeast
MW	16.9 kDa
Entrez	101693444
Expiration Date	6 months from date of receipt.