

Product Datasheet

Ferret IL-22 Recombinant Protein (orb1819626)

| | |
|-------------------------|---|
| Description | The Ferret IL-22 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th... |
| Target | IL-22 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | EAVPIGSHCR LANSDFEQPY ITNRTFMLAG EASLADNNTD VRLIGDTLFH GVNKEDHCYL MKEVLNFTLE EVLLPQSHRF QPVMQEVESF LTKLSNRLSQ CHIDSDDQHI QRNVQNLKDT VKKLGEDGEI KAIGELNLLF IALRDACI (148) |
| Purity | 98% |
| Source | Yeast |
| MW | 16.9 kDa |
| Entrez | 101693444 |
| Expiration Date | 6 months from date of receipt. |