

Product Datasheet

Human CTLA-4 Recombinant Protein (orb1819627)

Description	The Human CTLA-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CTLA-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	KAMHVAQPAV VCLASSRGIAS FVCEYASPGK ATEVRVTVLR QADSQVTEVC AATYMMGNEL TFLDDSICTG TSSGNQVNLT IQGLRAMDTG LYICKVELMY PPPYYLGIGN GTQIYVIDPE PCPDS (126)
Purity	98%
Source	Yeast
MW	13.5 kDa
Entrez	1493
Expiration Date	6 months from date of receipt.