

Product Datasheet

Bat IL-4 Recombinant Protein (orb1819630)

Description	The Jamaican fruit-eating Bat IL-4 yeast-derived recombinant protein is not tagged and is naturally ...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	RKPNLTLEET IKTLNILTAR KDPRCMELPV ADVLAAPKNT TERETFCAA TALRQLYRHH ECASRRFLSR LDRSLSLAG QNTCATDGVR KSALKDFLER LRRIMQEKYA RG (112)
Purity	98%
Source	Yeast
MW	12.8 kDa
Entrez	119037707
Expiration Date	6 months from date of receipt.