

## Product Datasheet

# Squirrel Monkey IFN beta Recombinant Protein (orb1819636)

<b>Description</b>	The Squirrel Monkey IFN beta yeast-derived recombinant protein is not tagged and is naturally endoto...
<b>Target</b>	IFN beta
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MSYNLLGFLQ RRSNFQSQKL LRQLNGRLEY CLKDRMDFEI PKEITQPQRL QKEDVALIVY EMLQNIFTIF NRDFSRTGWN KTIVENFLAN VYQHIEHLKT ILKEKWEKKE FTREKSMSSL HLKRYYGRIL HYLMAKEYSH CAWTIVSAEI LRNFYFIDRL TGYL (164)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	20.0 kDa
<b>Entrez</b>	<a href="#">101035681</a>
<b>Expiration Date</b>	6 months from date of receipt.