

Product Datasheet

Bat IFN Lambda 3 Recombinant Protein (orb1819637)

Description	The Common Vampire Bat IFN lambda 3 yeast-derived recombinant protein is not tagged and is naturally...
Target	IFN Lambda
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	PPLLTLSMP SLVSPSINTD LPVGCTLVLM LMTTVLTRTA AVPVPTPLSA LPGARGCVVA QFKFLAPQDMKAFRRAKDTLEELLKPNRSCSSRPFRTRDLRQLQVWERPV EAELALTLKVLGSIANSTLGDILDQPLHMLRYIHTKLOACVPAQATAGPRPRGHLH HWLHRLQEASKKESTGCLEASVTFNLFRLTHDLQCVASGGLCI (214)
Purity	98%
Source	Yeast
MW	23.4 kDa
Entrez	112320026
Expiration Date	6 months from date of receipt.