

Product Datasheet

Bovine Regakine-1 Recombinant Protein (orb1819639)

Description	The Bovine Regakine-1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	Regakine
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	EANEEPAGNM RVCCFSSVTR KIPLSLVKNY ERTGDKCPQE AVIFQTRSQR SICANPGQAW VQKYIEYLDQ MSK (73)
Purity	98%
Source	Yeast
MW	8.3 kDa
Entrez	504773
Expiration Date	6 months from date of receipt.