

Product Datasheet

Ferret IFN Beta Recombinant Protein (orb1819641)

Description	The Ferret IFN beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	IFN Beta
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AMNYNLLRFQLSSSSVEQCPELLQLNGTSKDCLKDRMNFKIPPEIQKSQKFQKED IVLVTLEMFQKTSDIFRRNLSSTGWNESIVENLLATLHWQKEHLEEILEDIMQEEN VTWDHRTLLHLKRYLRLIVRYLKAKEFSVCAWTIVQAEILRSFFFLDKLTDSLPN (166)
Purity	98%
Source	Yeast
MW	19.8 kDa
Entrez	101679100
Expiration Date	6 months from date of receipt.