

Product Datasheet

Human IL-36 alpha Recombinant Protein (orb1819644)

Description	The Human IL-36 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-36 alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	KIDTPQQGSI QDINHRVWVL QDQTLIAVPR KDRMSPVTIA LISCRHVETL EKDRGNPIYL GLNGLNLCLM CAKVGDQPTL QLKEKDIMDL YNQPEPVKSF LFYHSQSGRN STFESVAFPG WFIAVSSEGG CPLILTQELG KANTTDFGLT MLF (153)
Purity	98%
Source	Yeast
MW	17.1 kDa
Entrez	27179
Expiration Date	6 months from date of receipt.