

Product Datasheet

Ferret IFN Lambda 3 Recombinant Protein (orb1819645)

Description	The Ferret IFN lambda 3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-f...
Target	IFN Lambda
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VPGPKPLKVF PDAGGCHLAQ FQSLSPQELQ AFKKAKDTFE ESLSLKAWSC RPRLFRTWD LQQLQVGERP VALEAEVALT LKVLETMADR SLGSILDQPL HTLRHIQSEL QACVEAQPLA GPRPRGRLHH WLHRLHEAPK KEPLGCLENS VMFNLFRLT RDLKCVASGD LCA (173)
Purity	98%
Source	Yeast
MW	19.4 kDa
Entrez	101672127
Expiration Date	6 months from date of receipt.