

Product Datasheet

Mouse VEGF-C Recombinant Protein (orb1819647)

Description	The Mouse/Rat/Hamster VEGF-C yeast-derived recombinant protein is not tagged and is naturally endoto...
Target	VEGF-C
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AAHYNTEILK SIDNEWRKTQ CMPREVCIDV GKEFGAATNT FFKPPCVSVY RCGGCCNSEG LQCMNTSTGY LSKTLFEITV PLSQGPKPVT ISFANHTSCR CMSKLDVYRQ VHSIIRR (117)
Purity	98%
Source	Yeast
MW	13.1 kDa
Entrez	22341
Expiration Date	6 months from date of receipt.