

Product Datasheet

Ferret IL-1RA Recombinant Protein (orb1819649)

| | |
|-------------------------|--|
| Description | The Ferret IL-1 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturally... |
| Target | IL-1RA |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | HPSGKRPCR M QVFRWDINQ KTFYLRNNQL VAGYLQGSNT KLEEKLDVVP IEPHTVFLGI HGGKLCCLSCV KSGDETRLQL EAVNITDLSK SKDQDKRFTF ILSDSGPTTR FESAACPGWF LCTALEADRP VSFTNTPEEA RTVTKFYFQK D (151) |
| Purity | 98% |
| Source | Yeast |
| MW | 17.1 kDa |
| Entrez | 101678718 |
| Expiration Date | 6 months from date of receipt. |