

Product Datasheet

Ferret IL-1RA Recombinant Protein (orb1819650)

Description	The Ferret IL-1 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturally...
Target	IL-1RA
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HPSGKRPCR M QVFRWDINQ KTFYLRNNQL VAGYLQGSNT KLEEKLDVVP IEPHTVFLGI HGGKLCCLSCV KSGDETRLQL EAVNITDLSK SKDQDKRFTF ILSDSGPTTR FESAACPGWF LCTALEADRP VSFTNTPEEA RTVTKFYFQK D (151)
Purity	98%
Source	Yeast
MW	17.1 kDa
Entrez	101678718
Expiration Date	6 months from date of receipt.