

Product Datasheet

Ferret CCL8 Recombinant Protein (orb1819651)

Description	The Ferret CCL8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QPGSVSIPIT CCFNVVKKRV PLQKLESYTR VTNSQCPQEA VIFKTKASKE ICADPQLKWV QDNMKHLDQK SSTQKP (76)
Purity	98%
Source	Yeast
MW	8.6 kDa
Entrez	101686814
Expiration Date	6 months from date of receipt.