

Product Datasheet

Human IFN beta Recombinant Protein (orb1819661)

Description	The Human IFN beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	IFN beta
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MSYNLLGFLQ RSSNFQCQKL LWQLNGRLEY CLKDRMNFDI PEEIKLQQF QKEDAALTY EMLQNIFAIF RQDSSSTGWN ETIVENLLAN VYHQINHLKT VLEEKLEKED FTRGKLMSSL HLKRYYGRIL HYLKAKEYSH CAWTIVRVEI LRNFYFINRL TGYLRN (166)
Purity	98%
Source	Yeast
MW	20.0 kDa
Entrez	3456
Expiration Date	6 months from date of receipt.