

## Product Datasheet

### Ferret CCL16 Recombinant Protein (orb1819663)

<b>Description</b>	The Ferret CCL16 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
<b>Target</b>	CCL16
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	QSKSPASLNS SHTCCLKHHE KVMPRKMVVG YRRALNCYLP AIIFVTKKNR EICANPQNKW VQEYLRDPNI PLLPPRNLAQ ARSVRP (86)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	9.9 kDa
<b>Entrez</b>	<a href="#">101674669</a>
<b>Expiration Date</b>	6 months from date of receipt.