

## Product Datasheet

# Guinea Pig IFN alpha 17 Recombinant Protein (orb1819666)

<b>Description</b>	The Guinea Pig IFN alpha 17 yeast-derived recombinant protein is not tagged and is naturally endotox...
<b>Target</b>	IFN alpha 17
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	CTLPQTHSQG NTRALILLEQ MRRISPFSC L KHRQDFGLAQ LDFEGKQVQK AQALSILHEM TRQTFNLFIS EDSSAAWNKS LLHSFCTGLH QQLNDLEACL TQETGLEEPP LMHEDSRLVV RKYFHRIAVY LKEKSYSPCA WEVVRAQIVR SLSASEKLAR KMQD (164)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	18.9 kDa
<b>Entrez</b>	<a href="#">100717234</a>
<b>Expiration Date</b>	6 months from date of receipt.