

Product Datasheet

Swine Prolactin Recombinant Protein (orb1819668)

Description	The Swine Prolactin yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	Prolactin
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LPICPSGAVN CQVSLRDLFD RAVILSHYIH NLSSEMFNEF DKRYAQGRGF ITKAINSCHT SSLSTPEDKE QAQQIHHEVL LNLILRVLRS WNDPLYHLVT EVRGMQEAPD AILSRAIEIE EQNKRLLEGM EKIVGQVHPG IKENEVYSVW SGLPSLQMAP EDTRLFAFYN LLHCLRRDSH KIDNYLKLLK CRIIYDSNC (199)
Purity	98%
Source	Yeast
MW	23.0 kDa
Entrez	396965
Expiration Date	6 months from date of receipt.