

Product Datasheet

Mallard IL-17A Recombinant Protein (orb1819670)

Description	The Mallard IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	KVIRPGLKPE SLFKQADAGC LTQKDAKFPQ TVRVNISISN MNQDNKVTLD ISKRSLSPWD YRIDEDHNRF PQVIADAECR YSRCVNSAGQ LDHSLNSIPI NQEILVLRRE LKGCHHSYRL EKKIITVGCT CVTPLIQHQA (140)
Purity	98%
Source	Yeast
MW	15.9 kDa
Entrez	101795384
Expiration Date	6 months from date of receipt.