

## Product Datasheet

# Cynomolgus Monkey CTLA-4 Recombinant Protein (orb1819682)

<b>Description</b>	The Cynomolgus Monkey sCTLA-4 yeast-derived recombinant protein is not tagged and is naturally endot...
<b>Target</b>	CTLA-4
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	KAMHVAQPAV VLANSRGIAS FVCEYASPGK ATEVRVTVLR QADSQVTEVC AATYMMGNEL TFLDDSICTG TSSGNQVNLT IQGLRAMDTG LYICKVELMY PPPYMGIGN GTQIYVIDPE PCPDSD (126)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	13.5 kDa
<b>Entrez</b>	<a href="#">102115124</a>
<b>Expiration Date</b>	6 months from date of receipt.