

Product Datasheet

Ferret IFN alpha Recombinant Protein (orb1819683)

Description	The Ferret IFN alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IFN alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	CDLPQDHGLL NWRALMLLRQ MRRVSTFSCL RDRNDFAFPQ EVFDGKPLQK AQAISVLHVT NQKTFHLFCT EASPAPWNTT LLEQLCSGLS EQMSHLAACP LQEAGVGETP LVNGDSILRN YFQRISLYLQ EKQYSHCAWE MVRAEIMKVF SSSTTLQERL RGKKQDLVQH GDDSH (175)
Purity	98%
Source	Yeast
MW	20.0 kDa
Entrez	101686116
Expiration Date	6 months from date of receipt.