

Product Datasheet

Ferret IFN alpha Recombinant Protein (orb1819684)

| | |
|-------------------------|--|
| Description | The Ferret IFN alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free... |
| Target | IFN alpha |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | CDLPQDHGLL NWRALMLLRQ MRRVSTFSCL RDRNDFAFPQ EVFDGKPLQK AQAISVLHVT NQKTFHLFCT EASPAPWNTT LLEQLCSGLS EQMSHLAACP LQEAGVGETP LVNGDSILRN YFQRISLYLQ EKQYSHCAWE MVRAEIMKVF SSSTTLQERL RGKKQDLVQH GDDSH (175) |
| Purity | 98% |
| Source | Yeast |
| MW | 20.0 kDa |
| Entrez | 101686116 |
| Expiration Date | 6 months from date of receipt. |