

## Product Datasheet

# Ferret IFN alpha Recombinant Protein (orb1819685)

<b>Description</b>	The Ferret IFN alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
<b>Target</b>	IFN alpha
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	CDLPQDHGLL NWRALMLLRQ MRRVSTFSCL RDRNDFAFPQ EVFDGKPLQK AQAISVLHVT NQKTFHLFCT EASPAPWNTT LLEQLCSGLS EQMSHLAACP LQEAGVGETP LVNGDSILRN YFQRISLYLQ EKQYSHCAWE MVRAEIMKVF SSSTTLQERL RGKKQDLVQH GDDSH (175)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	20.0 kDa
<b>Entrez</b>	<a href="#">101686116</a>
<b>Expiration Date</b>	6 months from date of receipt.