

Product Datasheet

Bat CCL3 Recombinant Protein (orb1819697)

Description	The Black flying fox Bat CCL3 yeast-derived recombinant protein is not tagged and is naturally endot...
Target	CCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SYDENKTNVA STCCFTYTSR PIPRGRVVTY FKTSGQCSKH GVIFLTKKGQ HICADPHNAW VQEYISHLEN PQLVRWSSVP EGSKEQK (87)
Purity	98%
Source	Yeast
MW	9.9 kDa
Entrez	102894292
Expiration Date	6 months from date of receipt.